

# Mixers and Blenders

## WARING BLENDER WITH GLASS CONTAINER

#152-60 MODEL 700G, 115-VOLT

#152-60-1 MODEL 800G, 230-VOLT

The Model 700G/800G is a single-speed blender with a one-liter heat-resistant glass container with handle and two-piece vinyl and styrene lid. With a no-load test speed of 22,000 RPM, it is perfect for high-shear blending.

**Size:** 8" x 9" x 15.5" (20 x 23 x 39 cm)

**Weight:** 10 lb 12 oz (4.9 kg)

### Components:

#152-64 Container, 1 L, Glass



#152-60 Model 700G Waring Blender

**2-SPEED HAMILTON BEACH® BLENDER, 32 OZ (946 ML) STAINLESS STEEL CONTAINER, 115-VOLT, 60 HZ**

#152-04

**Size:** 13" x 10.8" x 8" (33 x 27 x 20 cm)

**Weight:** 9 lb 4 oz (4.2 kg)

*Did you know?*

OFITE has an engineering services department ready to help you with your next custom design project.

**WARING BLENDER WITH EXCITER GEAR, MAGNETIC PICK UP, AND STAINLESS STEEL CONTAINER, 2-SPEED**

#122-210 115-VOLT

#122-209 230-VOLT

**Size:** 13.5" x 9.5" x 8.25" (34 x 24 x 21 cm)

**Weight:** 12 lb 13 oz (5.8 kg)



#122-209 Waring Blender

### Component:

#122-208-1 Exciter Gear, Ten-Tooth

#122-209-1 Waring Blender (230-Volt Only)

#122-210-1 Waring Blender (115-Volt Only)

#122-210-2 Computer Plug

#122-210-3 DBA Connector, 9-Pin

#122-210-4 Connector Shield

**2-SPEED HAMILTON BEACH® BLENDER WITH ½ GALLON (1.9 L) STAINLESS STEEL CONTAINER, MODEL HBF400**

#152-05 115-VOLT

#152-05-1 230-VOLT

**Size:** 18" x 7" x 8" (45 x 18 x 20 cm)

**Weight:** 9 lb 7 oz (4.3 kg)



#152-05 2-Speed Hamilton Beach® Blender