

Mixers and Blenders

WARING BLENDER WITH GLASS CONTAINER

#152-60 MODEL 700G, 115-VOLT

#152-60-1 MODEL 800G, 230-VOLT

The Model 700G/800G is a single-speed blender with a one-liter heat-resistant glass container with handle and two-piece vinyl and styrene lid. With a no-load test speed of 22,000 RPM, it is perfect for high-shear blending.

Size: 8" x 9" x 15.5" (20 x 23 x 39 cm)

Weight: 10 lb 12 oz (4.9 kg)

Components:

#152-64 Container, 1 L, Glass



#152-60 Model 700G Waring Blender

2-SPEED HAMILTON BEACH® BLENDER, 32 OZ (946 ML) STAINLESS STEEL CONTAINER, 115-VOLT, 60 HZ

#152-04

Size: 13" x 10.8" x 8" (33 x 27 x 20 cm)

Weight: 9 lb 4 oz (4.2 kg)

Did you know?

OFITE has an engineering services department ready to help you with your next custom design project.

WARING BLENDER WITH EXCITER GEAR, MAGNETIC PICK UP, AND STAINLESS STEEL CONTAINER, 2-SPEED

#122-210 115-VOLT

#122-209 230-VOLT

Size: 13.5" x 9.5" x 8.25" (34 x 24 x 21 cm)

Weight: 12 lb 13 oz (5.8 kg)



#122-209 Waring Blender

Component:

#122-208-1 Exciter Gear, Ten-Tooth

#122-209-1 Waring Blender (230-Volt Only)

#122-210-1 Waring Blender (115-Volt Only)

#122-210-2 Computer Plug

#122-210-3 DBA Connector, 9-Pin

#122-210-4 Connector Shield

2-SPEED HAMILTON BEACH® BLENDER WITH ½ GALLON (1.9 L) STAINLESS STEEL CONTAINER, MODEL HBF400

#152-05 115-VOLT

#152-05-1 230-VOLT

Size: 18" x 7" x 8" (45 x 18 x 20 cm)

Weight: 9 lb 7 oz (4.3 kg)



#152-05 2-Speed Hamilton Beach® Blender