

HOT PLATE WITH THERMOSTAT

#168-01 115-VOLT

#168-01-1 230-VOLT

Size: 6" x 6" x 6" (15 x 15 x 15 cm)

Weight: 2 lb (907 g)



#168-01 Hot Plate with Thermostat

HOT PLATE

#168-03 115-VOLT

#168-03-1 230-VOLT

Size: 8" x 7.5" x 6" (20 x 19 x 15 cm)

Weight: 3 lb 5 oz (1.5 kg)



#168-03 Hot Plate

SPIN DRYER FOR CUTTINGS

#152-80

Size: 15" x 8" x 7" (38 x 20 x 18 cm)

Weight: 4 lb (1.8 kg)

The Cuttings Spin Dryer is designed for centrifugal quick-drying of cuttings samples without damaging them. Its compact design makes it ideal for remote locations. It is constructed of stainless steel and brass.

Components:

- #152-81 Inner Dryer Basket
- #152-82 Outer Cup, Stainless Steel
- #152-83 Bearing Post
- #152-84 Rubber Drive Clutch
- #152-85 Motor Shaft Collar
- #152-86 Bearing
- #152-87 Motor, 1/5 HP, 115-Volt
- #152-88 Bracket Set
- #152-89 Microswitch
- #170-29 Power Cord with Male Plug, 6'

ACCESSORIES FOR MIXERS, BLENDEES, AND STIRRERS

- #152-37 AC Power Cord, 3-Conductor
- #152-37-1 AC Power Cord, 3-Conductor, Beige
- #152-38 AC Power Cord, 3-Conductor, International
- #152-47 Agitator for Hamilton Beach® Mixers, Butterfly
- #152-41 Agitator Button for Hamilton Beach® Mixers, Upper
- #152-42 Agitator Button for Hamilton Beach® Mixers, Lower
- #152-01 Armature for Model 936 Hamilton Beach® Mixer, 115-Volt
- #152-02 Armature for Model 936 Hamilton Beach® Mixer, 230-Volt
- #122-200 Blending Assembly for Waring Blender, 1 qt
- #122-207 Blade for Blender, 1 qt
- #152-03 Brush and Spring Set for Hamilton Beach® Mixer, 115-Volt
- #152-03-1 Brush and Spring Set for Hamilton Beach® Mixers, 230-Volt
- #152-51 Case for Multi Mixer®, Sound Proof, Stainless Steel
- #152-40 Container for Mixers, 30 oz, Stainless Steel
- #152-39 Container with Handle for Mixers, Stainless Steel, 30 oz
- #152-64 Container Assembly with Lid for Model 700, 1 L, Glass
- #152-07 Container Guide for Model 936
- #152-09 Container Guide for Model 950
- #152-06 Container Rest for Model 936
- #163-09 Impeller for 2-Speed Laboratory Mixer
- #152-43 Impeller Blade for Multi-Mixer
- #152-11 Motor for Models 936 and 950 Hamilton Beach® Mixers, 115-Volt
- #152-12 Motor for Models 936 and 950 Hamilton Beach® Mixers, 230-Volt
- #152-35 Powerstat, 115-Volt
- #152-36 Powerstat, 230-Volt
- #152-44 Screw for Lower Impeller for Hamilton Beach® Mixer
- #152-41-1 Screw for Bottom Agitator for Multi-Mixer
- #152-08 Screw for Container Guide for Model 936
- #152-53 Spacer for One-Agitator Spindle
- #153-53-11 Stirring Bar, Magnetic, 1/2" x 5/16"
- #153-53-6 Stirring Bar, Magnetic, 7/8" x 5/16"
- #153-53-1 Stirring Bar, Magnetic, 1" x 5/16"
- #153-53-4 Stirring Bar, Magnetic, 1 1/2" x 5/16"
- #153-53-5 Stirring Bar, Magnetic, 2" x 5/16"
- #153-53-8 Stirring Bar, Spin Wedge, 1 3/4" x 1/2"
- #130-74 Transformer, 230 to 115 Volts, 1.10 Amps, 50 / 60 Hz



#152-80 Spin Dryer for Cuttings

Optional:

- #130-74 Transformer, 230 to 115 Volt, 1.10 Amps, 50 / 60 Hz