

THREE SPINDLE HAMILTON BEACH® MIXER WITH CONTAINERS, THREE SPEED

#152-20 115-VOLT, 60 HZ
 #152-30 230-VOLT, 50 HZ

Size: 15" x 15" x 24" (38 x 38 x 61 cm)
 Weight: 39 lb (17.7 kg)

Components:

#152-40 Container, Stainless Steel



#152-20 Three Spindle Hamilton Beach® Mixer

CUP MIXER

#163-10 12-VOLT
 #163-20 115-VOLT

Size: 9" x 5" x 5" (23 x 12 x 12 cm)
 Weight: 2 lb 14 oz (1.3 kg)

Components:

- #130-46 Male Plug Connector (12-Volt Only)
- #152-42 Agitator Button, Lower
- #152-87 Motor, 1/15 HP (115-Volt Only)
- #163-10-001 Motor, 1/35 HP (12-Volt Only)
- #163-10-1 Power Cord, 18-2 (12-Volt Only)
- #163-20-1 Splash Guard
- #163-20-2 Shaft, Stainless Steel
- #163-20-3 Power Cord with In-line, Single-Speed Switch, 18-3 SJT, 8' (115-Volt Only)
- #163-20-4 All Thread, 8-32 x 4 1/2"



#163-10 Cup Mixer

SINGLE SPINDLE HAMILTON BEACH® MIXER WITH CONTAINER, THREE SPEED#152-00

115-VOLT, 60 HZ
 #152-10 230-VOLT, 50 HZ

Size: 7.25" x 8.5" x 21.25" (18 x 22 x 54 cm)
 Weight: 13 lb (5.9 kg)

Components:

#152-40 Container, Stainless Steel

Optional:

#152-00-SP Spare Parts for One Year for #152-00



#152-00 Single Spindle Hamilton Beach® Mixer

Did you know?

We want to make ordering easy for you! You can email us at sales@ofite.com, call us at 1-877-TEST-MUD (1-877-837-8683) or 713-880-9885, fax us at 713-880-9886, or shop online at www.ofite.com.