

HTHP Filtration

BACK PRESSURE RECEIVER FOR CO₂, 15 ML, STAINLESS STEEL TUBE

#170-06

Components:

- #143-00 Regulator, CONCOA/AIRCO
- #143-01 Gauge, 200 PSI, 1/8" Bottom Connection
- #143-02-10 CO₂ Puncture Head Assembly
- #143-03 Barrel for CO₂ Bulb
- #143-06 Safety Bleeder Valve, 1/4" NPT
- #144-11 Street Ell, 1/8"
- #170-07 O-ring for Receiver
- #170-28 Receiver Body, 15 mL, Stainless Steel
- #170-32 Needle Valve, Male, 1/8" x 1/8" NPT
- #171-22 Retainer Pin



#170-06 Back Pressure Receiver

BACK PRESSURE RECEIVER FOR N₂, 15 ML, STAINLESS STEEL TUBE

#170-06-1


Components:

- #144-11 Street Ell, 1/8"
- #144-15 Bushing, 1/4" NPT Male to 1/8" NPT Female, Plated Brass
- #170-07 O-ring
- #170-28 Receiver Body
- #170-32 Needle Valve, Male, 1/8" x 1/8" NPT
- #171-22 Retainer Pin

BACK PRESSURE RECEIVER, 100 ML #171-10

Components:

- #170-32 Needle Valve, Male, 1/8" x 1/8" NPT
- #171-11 O-ring for Receiver Body
- #171-12 Receiver Body
- #171-22 Retainer Pin



#171-10 Back Pressure Receiver, 100 mL

