

HTHP Filtration

HTHP SINGLE-END CELL ASSEMBLY FOR CEMENT TESTING, 500 ML, 10" (25.4 CM), 316 STAINLESS STEEL, 2,000 PSI #171-18

Size: 12" x 10" x 5" (31 x 25 x 13 cm)
Weight: 16 lb 8 oz (7.5 kg)

Components:

- #170-13 O-ring for Cell, Buna N
- #170-16 Valve Stem
- #170-17 O-ring for Valve Stem, Viton®
- #170-18 Screen, 325-Mesh with 60-Mesh Backup, Detachable
- #170-24 End Cap with Removable Screens, 2,000 PSI
- #170-26 Cap Locking Screw, Stainless Steel
- #170-27 Allen Wrench for Cap Locking Screw, 5/16"
- #171-17 Cell Body, 10", 500 mL

HTHP SINGLE-END CELL ASSEMBLY FOR FILTER PAPER, 500 ML, 10" (25.4 CM), 303 STAINLESS STEEL, 2,000 PSI #171-20

Size: 12" x 10" x 5" (31 x 25 x 13 cm)
Weight: 16 lb 8 oz (7.5 kg)

Components:

- #170-13 O-ring for Cell, Buna N
- #170-16 Valve Stem
- #170-17 O-ring for Valve Stem, Viton®
- #170-26 Cap Locking Screw, Stainless Steel
- #170-27 Allen Wrench for Cap Locking Screw, 5/16"
- #171-17 Cell Body, 10", 500 mL
- #171-21 Cell Cap with 60-Mesh Screen, 2,000 PSI




#171-20 Single-End Cell Assembly, 500 mL

HTHP DOUBLE-END CELL ASSEMBLY FOR CERAMIC DISKS, 500 ML, 10" (25.4 CM), 316 STAINLESS STEEL, 2,000 PSI #171-29

Size: 12" x 10" x 5" (31 x 25 x 13 cm)
Weight: 16 lb 8 oz (7.5 kg)

Components:

- #170-13 O-ring for Cell, Buna N
- #170-16 Valve Stem
- #170-17 O-ring for Valve Stem, Viton®
- #170-26 Cap Locking Screw, Stainless Steel
- #170-27 Allen Wrench for Cap Locking Screw, 5/16"
- #170-69 Scribed Cell Cap for Ceramic Disk
- #170-72 Spacer for Filter Paper, 1/4", Stainless Steel
- #170-77 O-ring for Stainless Steel Spacer
- #171-02 Cell Body for Ceramic Disks, Double-End, 2,000 PSI, 500 mL
- #171-21 Cell Cap with 60-Mesh Screen, Inlet, 2,000 PSI



#171-29 Cell Assembly for Ceramic Disks, 500 mL

HTHP DOUBLE-END CELL ASSEMBLY FOR CEMENT TESTING, 500 ML, 10" (25.4 CM), 303 STAINLESS STEEL, 2,000 PSI #171-19

Size: 12" x 10" x 5" (31 x 25 x 13 cm)
Weight: 16 lb 8 oz (7.5 kg)

Components:

- #170-13 O-ring for Cell, Buna N
- #170-16 Valve Stem
- #170-17 O-ring for Valve Stem, Viton®
- #170-18 Screen, 325-Mesh with 60-Mesh Backup, Detachable
- #170-24 End Cap for Detachable Screen, 2,000 PSI
- #170-26 Cap Locking Screw, Stainless Steel
- #170-27 Allen Wrench for Cap Locking Screw, 5/16"



#171-19 Double-End Cell Assembly, 500 mL