


BACK PRESSURE RECEIVER FOR MODEL MB FILTER PRESSES

#171-58

Components:

- #142-37 Regulator, Victor
- #142-58 O-ring for Cell Coupling
- #142-61 Gauge, 200 PSI, 1/4" Bottom, 2" Face
- #143-02-10 CO₂ Puncture Head Assembly
- #143-03 Barrel for CO₂ Cartridge
- #143-06 Safety Bleeder Valve, 1/4" NPT
- #143-09 Relief Valve, 200 PSI (1,379 kPa), 1/4" NPT
- #170-07 O-ring for Receiver
- #170-32 Needle Valve, Male, Outlet, 1/8" x 1/8" NPT
- #171-66 Receiver Tube, Stainless Steel
- #171-75 Receiver Body




#171-58 Back Pressure Receiver

DUAL NITROGEN MANIFOLD, 1,350 PSI AND 750 PSI

#171-24

Components:

- #170-20 Manifold Block
- #170-32 Needle Valve, Male, 1/8" x 1/8" NPT
- #171-22 Retainer Pin
- #171-24-001 Modified Regulator, Low-Pressure Side
- #171-24-002 Modified Regulator, High-Pressure Side
- #171-24-1 Nut, L.H., Regulator Inlet CGA-580
- #171-24-2 Nipple with Filter for Regulator Inlet 15-3SF
- #171-24-3 Union Elbow, Female, 1/4" Flare x 1/8" FNPT
- #171-24-4 Pipe Plug, 1/4" NPT, 316 Stainless Steel
- #171-24-5 Street Tee, 1/4" NPT, 316 Stainless Steel
- #171-25-1 Relief Valve, 750 PSI (5,171 kPa)
- #171-25-2 Relief Valve, 1,350 PSI (9,308 kPa)
- #171-26 Hose, 3,000 PSI, 3/16" x 3'
- #171-28 Dual Manifold Body
- #171-38 Gauge, 1,000 PSI, 2 1/2" Face, 1/4" Bottom Connection
- #171-40 Gauge, 1,500 PSI, 2 1/2" Face, 1/4" Bottom Connection
- #171-42 Gauge, 3,000 PSI, 2 1/2" Face, 1/4" Bottom Connection
- #171-90-06 Reducing Bushing, 316 Stainless Steel
- #171-90-07 Hex Nipple, 1/4" NPT, 316 Stainless Steel
- #171-90-13 Adapter, 1/4" Flare x 1/4" Male NPT



#171-24 Dual Nitrogen Manifold

HIGH-PRESSURE NITROGEN PRESSURING ASSEMBLY WITH TANK

#171-31

Components:

- #143-06 Safety Bleeder Valve
- #144-17 Adapter, 1/4" FNPT x 1/4" MNPT, 316 Stainless Steel
- #170-20 Manifold Block
- #170-32 Needle Valve, Male, 1/8" x 1/8" NPT
- #170-37 Nitrogen Cylinder, Right-Hand Thread, 21" x 7"
- #171-22 Retainer Pin
- #171-24-1 Nut for Inlet, Right-Hand Thread, CGA-580, Chrome
- #171-24-2 Nipple with Filter for Regulator Inlet
- #171-24-3 Union Elbow, Female, 1/4" Flare x 1/8" FNPT, Chrome Plated
- #171-26-1 Hose, 5000#, 3/16" x 2'
- #171-38 Gauge, 1000 PSI, 2 1/2", 1/4" NPT Bottom
- #171-42 Gauge, 3000 PSI, 2 1/2", 1/4" NPT Bottom
- #171-53 Regulator, High-Pressure, Victor,
- #171-90-06 Reducing Bushing, 1/4" MNPT x 1/8" FNPT, 316 Stainless Steel



#171-31 High-Pressure Nitrogen Assembly