

UNIVERSAL HEAT CUP FOR MODEL 900 VISCOMETER

#130-76-10 115-VOLT

#130-76-10-1 230-VOLT

Size: 4.5" x 4.5" x 4" (11 x 11 x 10 cm)

Weight: 3 lb (1.36 kg)

Components:

#130-31 Thermostat, 50° - 300°F
 #130-38-10 Locator Pin
 #130-38-2 Red Lens for Lamp
 #130-38-3 Lamp
 #130-38-7 Cup, Stainless Steel
 #130-76-10-004 Sheath
 #130-76-10-005 Bottom
 #130-76-10-2 Heater, 150-Watt, 115 / 230-Volt
 #130-76-10-3 Sheath Cover
 #130-76-10-4 Power Cord
 #130-76-10-5 O-ring, Viton® 75
 #130-76-10-6 Strain Relief
 #130-76-10-7 Locknut
 #130-76-10-8 O-ring for Strain Relief
 #135-43 Base for Indicator Lamp
 #171-32 Midget Knob



#130-76-10 Universal Heat Cup

DOUBLE-WALLED CIRCULATING CUP #130-15

Size: 4.5" x 4" x 6" (11 x 10 x 15 cm)

Weight: 2 lb 3 oz (1 kg)

Components:

#130-15-02 O-ring, #040, N304
 #130-15-04 O-ring, #043, N304
 #130-38-10 Locator Pin
 #130-38-7 Cup, Stainless Steel
 #130-76-10-52 Screw, Socket Button Head, 6-32 x ¼", Stainless Steel



#130-15 Double Walled Circulating Cup

REFRIGERATED AND HEATED BATH/CIRCULATOR, CARON MODEL 2050W, 5.2 L

#152-55 115-VOLT, 60 HZ

#152-55-1 230-VOLT, 50 HZ

Size: 13.5" x 13.5" x 14" (34 x 34 x 36 cm)

Weight: 49 lb (22.2 kg)



#152-55 CARON Model 2050W Bath/Circulator

Did you know?

OFITE offers the most competitive labor rates in the industry!

Optional:

#152-55-2 Hoses and Fittings Set
 #152-55-3 Case
 #152-55-4 Ethylene Glycol, 1 Liter