

# Low-Pressure Filtration


## CO<sub>2</sub> PRESSURE ASSEMBLY WITH TOP CAP #142-00

**Size:** 9.5" x 2.5" x 5" (24 x 6 x 13 cm)  
**Weight:** 2 lb 9 oz (1.2 kg)

### Components:

- #141-02 Top Cap for Cell
- #141-05 Gasket, 3/32", Neoprene
- #141-22 Filter, Felt
- #143-00 Regulator, CONCOA
- #143-01 Gauge, 200 PSI, 1/8" Bottom Connection
- #143-02-10 CO<sub>2</sub> Puncture Head Assembly
- #143-03 Barrel for CO<sub>2</sub> Cartridge, Polished Chrome
- #143-06 Safety Bleeder Valve, 1/4" NPT

**Note:** CO<sub>2</sub> bulbs must be ordered separately.



#142-00 CO<sub>2</sub> Pressure Assembly with Top Cap



#142-10 CO<sub>2</sub> Pressure Assembly

## OFITE CO<sub>2</sub> PUNCTURE HEAD ASSEMBLY #143-02-10

### Components:

- #143-02-11 Puncture Pin Holder Assembly, OFITE Design
- #143-02-12 Puncturing Pin
- #143-02-13 O-ring for CO<sub>2</sub> Bulb
- #143-02-14 O-ring for Puncturing Pin Holder Assembly

**Note:** Must use #143-03 barrel with this puncture assembly.



#143-02-10 CO<sub>2</sub> Puncturing Head Assembly

## LOW-PRESSURE NITROGEN PRESSURING ASSEMBLY WITH TOP CAP AND TANK #171-31-1

### Components:

- #141-02 Top Cap for API Filter Press
- #141-14 Air Hose, Low Pressure, 3'
- #141-19 Air Hose Adapter
- #143-06 Safety Bleeder Valve
- #170-36 Regulator Assembly for Nitrogen, Victor, 200 and 3,000 PSI Gauges
- #170-37 Nitrogen Cylinder, Right-Hand Thread, 21" x 7"



#171-31-1 Low-Pressuring Nitrogen Assembly